

EXOSC4-exosome component 4 Gene

Reference: PTXBC002777

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RRP41; RRP41A; Rrp41p; SKI6; Ski6p; hRrp41p; p12A; exosome complex component RRP41; exosome complex exonuclease RRP41; exosome component Rrp41; ribosomal RNA-processing protein 41; exosome component 4		

For research use only.

GENE DESCRIPTION	exosome component 4
GENE ID	54512
ORIGIN SPECIES	Human
ACCESSIONS	BC002777
VECTOR	pDONR223
NCBI SYMBOL	EXOSC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggctggagctctgtcggaccagggtaccgggtggacggggcgcgcggggagctgcgaagatccaggcggatggcgtgttcgcgaggctgacggctcggcctacattgagcaggcaacacacaaaggcactggctgtgtctacggccgcagagatc
cgggctcccggctcgagccctgccggacaggccctagtgaaactgtcaatatagttcagcagcttcagcacaggtgagcgaagcagcggccacatggggaccgttaagtctgtgagatgggctgcagctcggccagacttcgaagcagccatcctcacacagctgcacca
cgctcccagattgatactatgtcagggtctacaggcagatggggacctatgcagctgtgtgaatgcagccagctggcagtgctggatccgggataccatgagagacttgtgtgtcgtcctcagctcgtcgtggacggcagacggccctggcgacctagccatgtggagg
aagcagctgtggccccagctggccctggcctgtgccagctcaggacagattgcgctgtgagatgatcccggctgcacgaggaccactggagcgggtgtggaggctgtgccaggctgccgagatgtgcacaccttagatcagtggtccggcagcatgtc
gtgaggcctctatctgtggggactga