



Reference: PTXBC002780

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FMO2; dimethylaniline monooxygenase [N-oxide-forming] 4; FMO 4; dimethylaniline oxidase 4; hepatic flavin-containing monooxygenase 4; flavin containing monooxygenase 4		

For research use only.

GENE DESCRIPTION	flavin containing monooxygenase 4
GENE ID	2329
ORIGIN SPECIES	Human
ACCESSIONS	BC002780
VECTOR	pDONR223
NCBI SYMBOL	FMO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]