



Reference: PTXBC002786

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEP27; SDR25C1; dehydrogenase/reductase SDR family member 2, mitochondrial; dehydrogenase/reductase (SDR family) member 2; dehydrogenase/reductase member 2; dicarbonyl reductase HEP27; protein D; short chain dehydrogenase/reductase family 25C member 1; sh		

For research use only.

GENE DESCRIPTION	dehydrogenase/reductase (SDR family) member 2
GENE ID	10202
ORIGIN SPECIES	Human
ACCESSIONS	BC002786
VECTOR	pDONR223
NCBI SYMBOL	DHRS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1