

RNF181-ring finger protein 181 Gene

Reference: PTXBC002803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF181; HSPC238; ring finger protein 181		
For research use only.			

GENE DESCRIPTION	ring finger protein 181
GENE ID	51255
ORIGIN SPECIES	Human
ACCESSIONS	BC002803
VECTOR	pDONR223
NCBI SYMBOL	RNF181
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcctatttcgatgaacacgactgcgagccgtcgaccctgagcaggagacgcgaaccaatgctgctggagctcgaaggtcacttttcaataggatggactttgaagacttgggttggttagatattgggaccaccacgtgcctccaccagctgccaagactgtggttgaga
acctccccaggacagtcacagagctctcagctcagctcaagtgccccgtgtgtcttttggaaattgaggaggagagactgccattgagatgccttgccatcaccttttccattccagctgcattctgccctggctaaagacaaattcctgtccctgtgcccgtatgagctgccactg
atgacgacacttatgaggagcacagacagagataaggctcgaacacagcagcagcaacacgactggagaacctcatggagccatgtacacgtga