

NOL12-nucleolar protein 12 Gene

Reference: PTXBC002808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Nop25; dJ37E16.7; nucleolar protein 12		
For research use only.			

GENE DESCRIPTION	nucleolar protein 12
GENE ID	79159
ORIGIN SPECIES	Human
ACCESSIONS	BC002808
VECTOR	pDONR223
NCBI SYMBOL	NOL12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcccaacaagaagaagcagatgggtgacggcgccggcgaggtctgttcttagcttcgacgaggagaagaggcgggagtagctgacaggctccacaagcggaaggtcgagcgaagaaggcagccattgaggagattaagcagcggtgaaaggagcagag
gaagcttcgggaggagcggccaccaggaatactgaagatgctggcagagagagaagggctctggaggaggcagatgagctggaccgggttggtgacagcaagacggagtcgggtgcatgaccacccaaccacacagtcacgtgaccacatcagtgacctgctcg
ggggcccggtgctcgggtgacccacctgagggaggggctggagacaggtctgaggaggaggcgtcatccacgggagaaccaaccaaaagccttgcccagggaagtcagagacccctgctctcagcgggatctctccctcacagcatcactacatgcacagcccaaaa
ggtcaaggagaaacatccccgacgggcccagactccaaaaagccccaaggggccctctgaccagcaaggcccagcgccgctctcacaggcaagcacggcacagcggggagtga