



Reference: PTXBC002810

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Rhith; Zfp13; zinc finger protein 205; transcriptional repressor Rhith; zinc finger protein 210		

For research use only.

GENE DESCRIPTION	zinc finger protein 205
GENE ID	7755
ORIGIN SPECIES	Human
ACCESSIONS	BC002810
VECTOR	pDONR223
NCBI SYMBOL	ZNF205
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr