

S100A2-S100 calcium binding protein A2 Gene

Reference: PTXBC002829

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAN19; S100L; protein S100-A2; S100 calcium binding protein A2		
For research use only.			

GENE DESCRIPTION	S100 calcium binding protein A2
GENE ID	6273
ORIGIN SPECIES	Human
ACCESSIONS	BC002829
VECTOR	pDONR223
NCBI SYMBOL	S100A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcagttctctggagcagcgctgctgtgctgctcactaccttcacaaagtactctgccaagagggcgacaagtcaagctgagtaagggggaatgaaggaaactctgcacaaggagctgccagctttgtgggggagaaagtggatgaggaggggctgaagaagctgatggcagcctggatgagaacagtgaccagcaggtggactccaggagtatgctgttttctggcactcatcactgtcatgtgcaatgactttccagggtgccagaccgacctga