

NRGN-neurogranin (protein kinase C substrate, RC3) Gene

Reference: PTXBC002835

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RC3; hng; calmodulin-binding protein; neurogranin (protein kinase C substrate, RC3); protein kinase C substrate		

For research use only.

GENE DESCRIPTION	neurogranin (protein kinase C substrate, RC3)
GENE ID	4900
ORIGIN SPECIES	Human
ACCESSIONS	BC002835
VECTOR	pDONR223
NCBI SYMBOL	NRGN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactgctgaccgagacgcctgctccaagccggacgacattctagacatcccgctggacgatcccgccgccaacgcggccgccaataccaggcgagtttctggggccacatggcgcggaagaagataagagcggagagcgcggccggaagggccgggcccc
tggggggcctggcgagctgggggtggccggggaggcgcgggcgcccgccagcgagactag