

SAR1B-SAR1 homolog B (*S. cerevisiae*) Gene

Reference: PTXBC002847

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GTP-binding protein SAR1b; ANDD; CMRD; GTBPB; SARA2; 2310075M17Rik; GTP-binding protein B; GTP-binding protein Sara; SAR1 homolog B; SAR1a gene homolog 2; secretion associated Ras related GTPase 1B		

For research use only.

GENE DESCRIPTION	SAR1 homolog B (<i>S. cerevisiae</i>)
GENE ID	51128
ORIGIN SPECIES	Human
ACCESSIONS	BC002847
VECTOR	pDONR223
NCBI SYMBOL	SAR1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctcatatttgattggatttacagtggttcagcagtggtgtctacagtttttaggattatataagaaaactggtaactgggtatttcttgattggataatgcagaaaaacaacattgctacacatgctaaaagatgacagacttggaacaatgtcccaattacatccactccgaagaac
tgaccattgctggcatgacgtttacaactttgatctgggtggacatgttcaagctcgaagagtggtgaaaaactaccttccctgctatcaatggcattgtatttctggtgattgtgcagaccagaaaggctgttagagtcaaaagaagaactgattcactaatgacagatgaaacattgctaat
gtgccttactgattcttgggaataagatcgcagaccctgaagccatcagtgaaagagaggttgcgagagatgtttgtttatatggtcagacaacaggaaagggaggtatctctgaaagaactgaatgcccgacccttagaagtgttcatgtgtgtgtcctcaaaagacaaggttacggag
aaggcttccgctggatggcagatcattgattaa