

COPZ1-coatomer protein complex, subunit zeta 1 Gene

Reference: PTXBC002849

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | CGI-120; COPZ; HSPC181; zeta-COP; zeta1-COP; coatomer subunit zeta-1; coatomer protein complex, subunit zeta; zeta-1 COP; zeta-1-coat protein; coatomer protein complex subunit zeta 1 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | coatomer protein complex, subunit zeta 1 |
| GENE ID | 22818 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC002849 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | COPZ1 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atggaggcgtgatttggaaacctccctgtatactgtcaagccatcctgattctggacaatgatggagatcgacttttccaagtactatgacgacacctccccagtgcaaggagcaaaaggcctttgagaagaacatttcaacaagaccatcgactgacagtgaattgccctctg
gaagcctgacagtggtatacaaaagcagtagatctctatctatgattggcagctcctatgaaaatgagctgatcctatggctgttctgaaactgtctctgactcattgagccagatgctgaggaaaaatgagaaaaagcagcactgctgagaacatggaggcctgttctggct
gtggatgaaattgtagatggagggtgactctagagagtgatccccagcaggtggtacaccgggtggcattaagggggtgaagatgccccctacggagcagaccgtgtctcaggtgctgagtcagccaagaacagatcaagtggtcactcctcgggtga