

ZNF434-zinc finger protein 434 Gene

Reference: PTXBC002859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF434; HCCS-5; zinc finger and SCAN domain-containing protein 32; human cervical cancer suppressor gene 5 protein; zinc finger protein 434; zinc finger and SCAN domain containing 32		

For research use only.

GENE DESCRIPTION	zinc finger protein 434
GENE ID	54925
ORIGIN SPECIES	Human
ACCESSIONS	BC002859
VECTOR	pDONR223
NCBI SYMBOL	ZNF434
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcagtgaagagtaccgagccacacccatctcaaacaggatccacacaggccagaaatcagccctcagggttaacagccctgactccgagccctccgtcagcgttcaggcagtttctaccaggaggttaactggcccatgaagcttttagcaaaactctgggaactctgttgcagtggtgagccgaagaccactcaaaagaggaaatctggagctgctgttttgagcagtttctgactatcttccagaggagatccagacctgggtgagggagcagcatccagaaaaacggcaggaagctgtgctctgtgttgaggtatgacagagagctcctggacaacaggttctagatctgagaaggacttgaagtactcatgaaggagatggccctttgggagcaaccagagaaatcactgagatcccaatggaaacaggaggttcagccagaggaaccgacttttaagggtacacagagctcacaccaaaagaccagggaacagtcagaagcctggctgctctcaggtcccaggaaactgcctcaaacagggtctccagcaggagacaggtgctgtgtgtctgacagctgggtccaggaccagccatgcgtgacaacagagctgtatccctctgtcagcaagaatggtgtgccaggccctgcacaaagggccctctacaggggtgcacccagaggaaggacagtcacgtctcgtggcaacaggtgccctgggctatga