

ATAD4-ATPase family, AAA domain containing 4 Gene

Reference: PTXBC002865

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATAD4; proline-rich protein 15-like protein; ATPase family, AAA domain containing 4; testicular tissue protein Li 157; proline rich 15 like		

For research use only.

GENE DESCRIPTION	ATPase family, AAA domain containing 4
GENE ID	79170
ORIGIN SPECIES	Human
ACCESSIONS	BC002865
VECTOR	pDONR223
NCBI SYMBOL	ATAD4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgactgaaattggttggaagctgacttctccgaaaaagaatccactcccaaaagtgtgtatgagatccctgacacatgcccacagagggagatgcagaacccccgaggcctgacgctggagggcccaacagcgactttaacaccgcctggagaagattgtgga
caagagcacaagggcaagcacgtcaaggtctcgaactcaggacgttcaaggagaagaagaagtgtgagccacgctggcagagaacccctacctctttgatgacgaggaaggacggtcatcaaatga