



Reference: PTXBC002876

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD40 repeat-containing protein SMU1; BWD; SMU-1; fSAP57; brain-enriched WD-repeat protein; functional spliceosome-associated protein 57; smu-1 suppressor of mec-8 and unc-52 homolog; smu-1 suppressor of mec-8 and unc-52 protein homolog; DNA replication re		

GENE DESCRIPTION	smu-1 suppressor of mec-8 and unc-52 homolog (C. elegans)
GENE ID	55234
ORIGIN SPECIES	Human
ACCESSIONS	BC002876
VECTOR	pDONR223
NCBI SYMBOL	SMU1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1