

MED27-mediator complex subunit 27 Gene

Reference: PTXBC002878

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRAP34; CRSP34; CRSP8; MED3; TRAP37; mediator of RNA polymerase II transcription subunit 27; CRSP complex subunit 8; cofactor required for Sp1 transcriptional activation, subunit 8, 34kDa; p37 TRAP/SMCC/PC2 subunit; transcriptional coactivator CRSP34; med		

For research use only.

GENE DESCRIPTION	mediator complex subunit 27
GENE ID	9442
ORIGIN SPECIES	Human
ACCESSIONS	BC002878
VECTOR	pDONR223
NCBI SYMBOL	MED27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcggaacaaggagacgctggaggccgggagaagcccttattgcgcactccaggacaactacattcggtcaaccgggacctcaatgagctggaacgtctgagcaactctggtaggcaagccatctgagaaccatctctcataacagtgggctgttaagcctggatcctgtgcag
gacaaaactcctctctatagtcactcctttaa