

NECAP1-NECAP endocytosis associated 1 Gene

Reference: PTXBC002888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIEE21; adaptin ear-binding coat-associated protein 1; NECAP endocytosis associated 1		

For research use only.

GENE DESCRIPTION	NECAP endocytosis associated 1
GENE ID	25977
ORIGIN SPECIES	Human
ACCESSIONS	BC002888
VECTOR	pDONR223
NCBI SYMBOL	NECAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgctcgtcctaagttggatctgggttcaaggaaggacaacacatcaagttgtgtatcgggaacattacaacaagaaggaggtgcttctaagccaggactgcaaggggtgggtctgagcttactccaccccgccaggaggcaagtcactattccccaccatctctcagttgccatcagcaatcatgtcaccaccacccattccgaaatcaaccatggagcagtgatgcagatacttttagatttggattctctgctcctgtcacgacaccagcaccaactccagtttctgtaagcaatgacttggggagacttcagcactgcctccagctctgttccaaacagacacacagccatccaactgggtccagttctga