

TMEM160-transmembrane protein 160 Gene

Reference: PTXBC002907

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 160		
For research use only.			

GENE DESCRIPTION	transmembrane protein 160
GENE ID	54958
ORIGIN SPECIES	Human
ACCESSIONS	BC002907
VECTOR	pDONR223
NCBI SYMBOL	TMEM160
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaggcggctgggtgggctggggccgctcgtccgtcttcgttcgttcggaggtcgtactgcccgcctcagcggccccggagcggggcgccccgggggtcttcgccccggccacgggtcccgcgccgggggttcgccccggccaggtcccgagctggatcgtgcgg
acgcctggctcctcgaagagcgcagagacagccttcctcctggttcgcaatggcctcctggcatcgggcatcggggcctccttcacgagatgacatgggtcgggaagcagcatalggcttctcctgctggcgcgctgtgctgtgtgggcagcgcctcgtacgcgt
gagcctggcgcgctgcgaggacccatgcagctgacgtggggggcgccgctggggcgggcgccgtgctggccgcccagcctgcttggcgctgcggcgtgggcctctacatggggcagctggagctggacgtggagctgggtgccgaggacgacgggacggcctccgcg
gaaggccctgatgaggcgggtcggccgccacccgagtga