



Reference: PTXBC002927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSCARG; SDR48A1; nmrA-like family domain-containing protein 1; NmrA-like family domain containing 1; short chain dehydrogenase/reductase family 48A, member 1; NmrA like redox sensor 1		

For research use only.

GENE DESCRIPTION	NmrA-like family domain containing 1
GENE ID	57407
ORIGIN SPECIES	Human
ACCESSIONS	BC002927
VECTOR	pDONR223
NCBI SYMBOL	NMRAL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]