

B9D1-B9 protein domain 1 Gene

Reference: PTXBC002944

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPPB9; JBTS27; MKS9; MKSR1; B9 domain-containing protein 1; B9 protein domain 1; MKS1-related protein 1; endothelial precursor protein B9; B9 domain containing 1		

For research use only.

GENE DESCRIPTION	B9 protein domain 1
GENE ID	27077
ORIGIN SPECIES	Human
ACCESSIONS	BC002944
VECTOR	pDONR223
NCBI SYMBOL	B9D1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaccgcgagtcctagctcttctactcatggtcaacgggcaggtggagagcgcccggttcagagtagatgacgtctactgcaagtagctgttggtagcggcaggactgggccccacagcgggtctggaggaggggatctcacagatcacatccaagagccaagatgtg
cggcaagcactggtgtggaacttccccatgatgtcactttaaaagcaccacccctacggctggccacagatcgtctcagcgtgtatggaccagatgtgttcgggaacgatgtgttcgagggctatggggccgtgcacgtgcccctctcacctggccggcacaaaaggacatcccca
tggtttgccagaatctacgtctaaactgcagaagtttacaagtttggcctggcgcctcctcagatctgcaagctgctccaccactgaggacaaatag