



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD79B-CD79b molecule, immunoglobulin-associated beta Gene

Reference: PTXBC002975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD79b molecule; CD79b molecule, immunoglobulin-associated beta; CD79b antigen (immunoglobulin-associated beta); AGM6; B29; IGB; B-cell antigen receptor complex-associated protein beta chain; B-cell-specific glycoprotein B29; Ig-beta; immunoglobulin-associ		

For research use only.

GENE DESCRIPTION	CD79b molecule, immunoglobulin-associated beta
GENE ID	974
ORIGIN SPECIES	Human
ACCESSIONS	BC002975
VECTOR	pDONR223
NCBI SYMBOL	CD79B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagcggctggcgtgtgtctctgtgccagccactggatgggtggcgtgtgtctgtctgtctcagctgagccagttaccagcagccagatcggaggaccggtaccgggaalcccaaggtagtctgtttctgcggaatctggcagagccacgtttcatagccaggaaacggggcttcacg
gtgaaatgcatctgtcatgaacagcgctctccggcaatgtgagctgtgctctggaagcaggagatggagcagagaatccccagcagctgaagctgggaaaaagggccgcatgggaagatccagaaacgaatctctgcaccctcacctccaaggcatccggtttgaggacaatgatgcattc
cttcttcggagcagaaggtgcaacaacacactcggaggtgtgtctacacggctctggccagcagagctgcgagtcattgggattcagatcagcttcggcagatgctgaagcagctgaagcaggaacatcgtgaagatgggtatcatcatgatccagacgtctgtcatctcttcattcatctgtcctatctctctgtg

More info: Please contact us for any question: contact@proteogenix.fr