



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SPSB2-splA/ryanodine receptor domain and SOCS box containing 2 Gene

Reference: PTXBC002983

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRCC9; SSB2; SPRY domain-containing SOCS box protein 2; SPRY domain-containing SOCS box protein SSB-2; gene-rich cluster protein C9; splA/ryanodine receptor domain and SOCS box containing 2		

For research use only.

GENE DESCRIPTION	splA/ryanodine receptor domain and SOCS box containing 2
GENE ID	84727
ORIGIN SPECIES	Human
ACCESSIONS	BC002983
VECTOR	pDONR223
NCBI SYMBOL	SPSB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcagacagctctggcagggggagcagcagcagaccccacgccacagggccctctacacctgacctctctgccgaggccttggaagagctctgtctgtcacccccctctgacctgggggccccagcgcgcacagcttggaaccccaaaagactgttcagagaacatcagagt
caagggaaggagggtgttacttggagcggcgcccgctggcccgagagcactgatggggccgggggtaaggagggtctattcaaggggctgcacgcctgggagatcagctgccccctagagcagagggggcagcagatccgtgtgtggggctgtggccacggccctcgcccgctgcaga
ctgaccactacgcggcgctgtctggcagcaacagcagtgctgtggggctgggacatcgggcgaggagcctgtgatacatcagacgaaggggccgcgagcccccagatccagcgggaactcaggctgtagcagctgtgagggtgccagagacatgctgggtgtgttgacatggag
ggagggaactctggcctgctattggggccgacaccttggggccgacatctcggagggcagggcgactgaagggcagggccacctctcggcagtaagcctgtgtggggcagctggcagctccgcactcctggcgaaagagagcgggacccacctcccttctgacctgagccgcctg
tgtctgcgcacaacctctggcagataccgcgtcgcgcagctgtctgtccctgcccttccccctctccatgaagcctgacctctgacacatgta

More info: Please contact us for any question: contact@proteogenix.fr