

VAPA-VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa Gene

Reference: PTXBC002992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VAP-33; VAP-A; VAP33; hVAP-33; vesicle-associated membrane protein-associated protein A; 33 kDa VAMP-associated protein; VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa; VAMP-A; VAMP associated protein A		

For research use only.

GENE DESCRIPTION	VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa
GENE ID	9218
ORIGIN SPECIES	Human
ACCESSIONS	BC002992
VECTOR	pDONR223
NCBI SYMBOL	VAPA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaagcagcagcagatcctggtcctcgatccgccacagacctcaaatcaaggcccttcacagatgtagtactacaatcttaattgcgaatccatcggatagaaaagtgtgttcaaaagtgaagactacagcacctcgccggtactgtgtgagcccaacagtgaattatt
gacccagggtcaactgtgactgtttcagtaagtctacagcccttgactatgatccgaatgaaaagagtaaacacaaagtgttggtacagacaattttgctccacaaacactcagataggaaagctgtgtggaagaggcaaaacctgatgaattaatggatccaaatgagatgcgtatttg
aaatgcccaatgaaaatgataaattgaatgatatggaacctagcaaaagctgttccactgaatgcataagcaagatggacatgccaaaccacacagtggttacttaataccgaacaaaggaactaatggaagagtgtaaaagactcagggaagaatgatgaagctatcagaag
aaaatggcacctgagagatgaagggttaaggctcagaaggttagcacattcggataaacctggatcaacctcaactgcaccttcagagataatgtaccagtcctctcttcactcttctgttaattgcagccatttcattggattctttaggaaattcatctttag