

MRT04-mRNA turnover 4 homolog (*S. cerevisiae*) Gene

Reference: PTXBC003013

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ribosome assembly factor MRT04; MRT04 ribosome maturation factor; C1orf33; MRT4; dJ657E11.4; mRNA turnover protein 4 homolog; 60S acidic ribosomal protein PO; MRT4, mRNA turnover 4, homolog; mRNA turnover 4 homolog; MRT4 homolog, ribosome maturation facto		

For research use only.

GENE DESCRIPTION	mRNA turnover 4 homolog (<i>S. cerevisiae</i>)
GENE ID	51154
ORIGIN SPECIES	Human
ACCESSIONS	BC003013
VECTOR	pDONR223
NCBI SYMBOL	MRT04
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaaatccaagcgcacagaagaagtctcctaaccaaaactgccaaagagcttggaaatgaacaaaacctgatagaagagcttggaaatgtgtggacactacaagtacatttcatctctgtgtggccaacatgaggaaacagcaagctgaaggacatccggaacgcctgga
agcacagccggtatgttcttggcaaaaacaggtgatgatgtgtgacctgggtcgagagccatctgatgaatacaaaacacacacagctgacacaggtcagcaaaaggttgagggtgaggtgggtctctctgtttaccaaccgcacaaaggaaggtgaatgagtggttcacgaaatcacaca
gaaatggactacgcccagctggtaacaaagcagcttcaactgtgagcctggatccagggccctggagcagttccccactccatggagccacagctcaggcagctgggcctgccaccgccctcaagagaggtgtgtgactctgtctgtgactacgaggtgtgcaaggaggcg
atgtgtgacccagagcaggtcgcgtcgaagcttttgggtatgatggctgaattcaaggtgacatcaatacatgtgggattcacagtcgggaaggttcagcagatgggagacgacttgcagagagcgcactctgagtcacagaagagtcagactcagaagatgatgactg
a