

ARF5-ADP-ribosylation factor 5 Gene

Reference: PTXBC003043

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADP-ribosylation factor 5; ADP ribosylation factor 5		
For research use only.			

GENE DESCRIPTION	ADP-ribosylation factor 5
GENE ID	381
ORIGIN SPECIES	Human
ACCESSIONS	BC003043
VECTOR	pDONR223
NCBI SYMBOL	ARF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggectcaccgtgtccgcgctcttttcgcatcttcgggaagaagcagatcggattctcatgttgcttgatcggtggcaagaccacaatcgtacaaactgaagttggggagattgtcaccaccatccaacctaggttcaatgtagaacagtgaatataagaacatc
tgtttcacagctctgggacgtgggagccaggacaagattcggcctctgtggcggcactacttccagaacactcaggcctcatcttgggtggacagtaaatgaccgggagcgggtccaagaatctgtgatgaactccagaagatgctgcaggaggacgagctcgggatcagtgct
gctggattttgccaacaagcaggacatgcccacgcatgcccgtgagcagctgactgacaagctggggtacagcacttacgcagccgcacgtgggtatgtccaggccacctgtgccaccaaggcacaggtctgtacgatggctggactgctgccacgagctgtcaaacgct
aa