



Reference: PTXBC003045

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CED; CED1; SPG; WDR10p; WDR140; intraflagellar transport protein 122 homolog; WD repeat domain 10; WD repeat-containing protein 10; WD repeat-containing protein 140; intraflagellar transport 122 homolog; intraflagellar transport 122		

For research use only.

GENE DESCRIPTION	intraflagellar transport 122 homolog (Chlamydomonas)
GENE ID	55764
ORIGIN SPECIES	Human
ACCESSIONS	BC003045
VECTOR	pDONR223
NCBI SYMBOL	IFT122
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr