

ELMO1-engulfment and cell motility 1 Gene

Reference: PTXBC003051

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CED-12; CED12; ELMO-1; engulfment and cell motility protein 1; ced-12 homolog 1; engulfment and cell motility 1		

For research use only.

GENE DESCRIPTION	engulfment and cell motility 1
GENE ID	9844
ORIGIN SPECIES	Human
ACCESSIONS	BC003051
VECTOR	pDONR223
NCBI SYMBOL	ELMO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggtggtgaaggagcaggttatgagagcacttacaaccaagcctagctccctggaccagttcaagagcaaacctgcagacctgagctacactgagactcctgaaatccgagctccgagaggatgaaccaggaagattccagtcgccgcatgttggaactaaaggagaagatt
cagccagaaatcttagagctgatcaaacagcaacgctgaaccgcttgggaagggaacctgcttttaggaaactcaatgccggcgaggagcaagacaagtttgggtattgtcggcttccgcaaatcacaaagtcctgcatlacggagacttagaagagagtcctcagggaagtgcccc
acgattccttgaggacaaactgccgtggcagatataaagccgtggtgacgggaaaggactgccctcatatgaaagagaaggtgcccttaacaaacaaaggaggtgctgaaactcgtttcctcatctgtatgactcaaaactccaactgaactcatcgtctgacaaagcatgagt
actgtatctggacggatggactgaatgcctactcgggaaggacatgatgagcgacgtgacgcggaatgacctggacacctgctcagcatggaatcaagctccgctcctggacctggaaacatccagatccctgacgcacctccgcatccaagagcccgcaactatgac
ttcgtctatgactgaactga