



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CDK2-cyclin-dependent kinase 2 Gene

Reference: PTXBC003065

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	p33(CDK2); CDKN2; cyclin-dependent kinase 2; cdc2-related protein kinase; cell division protein kinase 2; p33 protein kinase; cyclin dependent kinase 2		

For research use only.

GENE DESCRIPTION	cyclin-dependent kinase 2
GENE ID	1017
ORIGIN SPECIES	Human
ACCESSIONS	BC003065
VECTOR	pDONR223
NCBI SYMBOL	CDK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtagaacctccaaaaggtgggaaaagatcggagagggcagctagcggagtggtgtgtacaaaagccagaaacaagttgacggggagaggtgtgtgtgcgttaagaaatccgcctggacactgagactgaggggtgtgccagctactgccatccgagagatctctgtcttaaggagcttaacca
tctaatattgtcaagctcgtggatgtcattcacacagaaaataaactctacctgtttttgaatttctgccaaagatctcaagaaatcatggaatgcctctgctctactggcattcctcttcccctcatcaagagctatctgttcacgtctgccagggcctgactttctgcattctatccgggtcct
caaccgacacttaaacctcagaatctgttatttaacacagagggggccatcaagctagcagacatcttggaactagcagagactctttggagtcctgtctgacttaccactgaggtgtgtgacccctgtgtaccgacacctggaatcctcttgggctcgaataattattccacagctgttgaca
cttggagcgtggctgcatctttgctgagctgacgtcggcggcctattctccctggagattctgagattgaccagctctccggactctggcagctctggcagccacagaggggtgtgtggcaggaggtacttacttactctgattacaacgaagtgtcccccaagtggcccccgaaga
ttttgaagtgtgtactccccctggtagaagtgagcgcagacttctatccgaatctctactacacaccttaaacagcgaattctggcgaagcagccctgcttaccctcttctcagagtgatccaggaacagcataccaccttctcagctctga