

CHCHD2-coiled-coil-helix-coiled-coil-helix domain containing 2 Gene

Reference: PTXBC003079

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf17; MNRR1; NS2TP; PARK22; coiled-coil-helix-coiled-coil-helix domain-containing protein 2; 16.7kD protein; HCV NS2 trans-regulated protein; aging-associated gene 10 protein; coiled-coil-helix-coiled-coil-helix domain-containing protein 2, mitochondri		

For research use only.

GENE DESCRIPTION	coiled-coil-helix-coiled-coil-helix domain containing 2
GENE ID	51142
ORIGIN SPECIES	Human
ACCESSIONS	BC003079
VECTOR	pDONR223
NCBI SYMBOL	CHCHD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgtggaagccgaagccgaactcccgcatggccctccggccagccgggcccctcagatgagagctgcaccagggcagcaccagtcgctcagccaccagcagggcaccctcatcgcagttggctctctgctgctgcgcccgccagccaggtctgatggccagat
ggcaaccactgcagctggcgtggctgtggctctgctgtggggcacacattgggtcacgccaacttactggggcttcagtgagggaagtaagctgagcctcgaggcctgacatcacttaccaggagcctcagggaaccagccagcacagcagcagcagccttgccctatgagatc
aaacagttctggagtggtccagaaccaggtgacatcaagctctgtgaggggttcaatgaggtgctgaacagtgccgacttcaaacggattggcctaa