



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

PLEKHJ1-pleckstrin homology domain containing, family J member 1 Gene

Reference: PTXBC003084

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GNRPX; pleckstrin homology domain-containing family J member 1; PH domain-containing family J member 1; guanine nucleotide releasing protein x; pleckstrin homology domain containing J1		

For research use only.

GENE DESCRIPTION	pleckstrin homology domain containing, family J member 1
GENE ID	55111
ORIGIN SPECIES	Human
ACCESSIONS	BC003084
VECTOR	pDONR223
NCBI SYMBOL	PLEKHJ1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcttacaacgagaagggagctcagggctctgtccggcgacccggcgagatggcgccgacgtcgggcatggggccccaagaagggcagcgtctcgaagcggcggtggtgaagctggtgtaattctcttctacttccggacagcagagccgacggccgtcggagccc
tgcctcgcgcgcgcagagctctcggcgagacggccgacgctctccatcagctcattgagacgctcgaggaagtatcatttgatgcagcagcaggagcagctgcaggagtgatggagcctctcgtcggcgagctacgagtcgagcggcgaggaagcctcatctctaca
ggaaacgaatccggaaggttgacggcgcaaggacccctcgaaacatttcgcatacggcagaggcagcttcagcttcgacgttcgagcttcgagctcgta

More info: Please contact us for any question: contact@proteogenix.fr