



Reference: PTXBC003089

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DXS8237E; GPATC9; GPATCH9; S1-1; TARPS; ZRANB5; RNA-binding protein 10; RNA-binding protein S1-1; g patch domain-containing protein 9; RNA binding motif protein 10		

For research use only.

GENE DESCRIPTION	RNA binding motif protein 10
GENE ID	8241
ORIGIN SPECIES	Human
ACCESSIONS	BC003089
VECTOR	pDONR223
NCBI SYMBOL	RBM10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]