



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PDLIM4-PDZ and LIM domain 4 Gene

Reference: PTXBC003096

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PDZ and LIM domain protein 4; LIM domain protein; LIM protein RIL; enigma homolog; reversion-induced LIM protein; PDZ and LIM domain 4		

For research use only.

GENE DESCRIPTION	PDZ and LIM domain 4
GENE ID	8572
ORIGIN SPECIES	Human
ACCESSIONS	BC003096
VECTOR	pDONR223
NCBI SYMBOL	PDLIM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacattcgtgacacctgcgcggcgctgcctctggggcttcgcctgtgctggcgccgggaactcagcgcgccctaccatctcagggtgcattcgtgcgacgaaggctcattggctgccctgtgccaggagacctgatccaggccatcaatggtgagagcacagagctcat
gacacacctggaggcacagaacccgacaaaggctgcacgatcactcacactgtctgtgagcaggcctgaggcgaggagctggcccatggtccctgatgacagaaggctcaggcacacaggatccacatgatctgagatccaggagggcagcccaacaaccaggcaggcg
gccctcaggcacccggactggccagaaatgacagacaaacgtggatctcatatggaacaacccctcgttttcagttctcacaatggcagcagcgaggcgccacctgtccaggccgaatgagcacctgcattgtctccacccccagcgtcagaccagccagagcgctccc
cgcgagcgccgggactcagatgcactggctgacgtgacatgtagtgcgtcgggagccggcgagccgctggcgccggcgagcccaagcagtcaggctctcttcgtactgtcaggcgatgctagaggccgcggaggcgccggggcaccatcgtcaaggcacgggacaagctcta
ccatcccgatgctcattgctgactgacgtgcgcctga