

LDOC1-leucine zipper, down-regulated in cancer 1 Gene

Reference: PTXBC003104

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein LDOC1; BCUR1; Mar7; Mart7; breast cancer, up-regulated 1; leucine zipper downregulated in cancer; leucine zipper protein down-regulated in cancer cells; leucine zipper down-regulated in cancer 1		

For research use only.

GENE DESCRIPTION	leucine zipper, down-regulated in cancer 1
GENE ID	23641
ORIGIN SPECIES	Human
ACCESSIONS	BC003104
VECTOR	pDONR223
NCBI SYMBOL	LDOC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggatgagttggtgctgctgcacgcgctcctgatgcggcaccgcgccctgagcatcgagaacagccagctcatggaacagctgcggctgctggtgtgcgagagggccagcctgctgcgccaggtacgtccgagctgcccgtgcccctcccgaaacgtttaatggcg
agagctccccgctccccgagtttatcgtgcagacggcgtcttacctgctcgtgaacgagaaccgattctgcaacgacgccatgaaggtggcattcctaatacgcctcctcaccggggaagccgaggagtggtggtgcctacatcgagatggatagccccatctaggtgattaccggg
ccttcctcgatgagatgaacagtccttggctgggatgacgacgaagacgacgacgacgaagaagaggaggatgattattag