



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FH-fumarate hydratase Gene

Reference: PTXBC003108

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

FMRD; HLRCC; LRCC; MCL; MCUL1; fumarate hydratase, mitochondrial; fumarase; fumarate hydratase

For research use only.

GENE DESCRIPTION	fumarate hydratase
GENE ID	2271
ORIGIN SPECIES	Human
ACCESSIONS	BC003108
VECTOR	pDONR223
NCBI SYMBOL	FH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1