

RPS3-ribosomal protein S3 Gene

Reference: PTXBC003137

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	40S ribosomal protein S3; IMR-90 ribosomal protein S3		
For research use only.			

GENE DESCRIPTION	ribosomal protein S3
GENE ID	6188
ORIGIN SPECIES	Human
ACCESSIONS	BC003137
VECTOR	pDONR223
NCBI SYMBOL	RPS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagtgcaaatccaagaagagggaagttgtcgtgatggcatcttcaagctgaatgaatttctactcgggagctggctgaagatggctactctggagttgaggtgcagttacacacacaggacagaatcattatcttagccaccagaacagaatgtttgttgaga
agggccggcggttcgggaactgactgctgtagttcagaagaggtttggcttccagagggcagtgtagagctttatgctgaaaagggtggccactagaggtctgtgtccattgcccaggcagagttctgtgttacaacctctaggaggctgtgtgcggaggcctgctatggtgtg
ctgcggttcatcatgagagtgggccaaggctgcgaggtgtgtgtgtctggaaactccgaggacagagggctaaatccatgaagttgtggatggcctgatgccacagcggagaccctgttaactactacgttgacactgctgtgcgccacgtgtgtctcagacaggggtgtgtgg
gcatcaaggtgaagatcatgctgcctgggacccaactggttaagattggcctaagaagccctgctgaccacgtgagcattgtgaacccaagatgagatactgccaccaccccatctcagaacagaagggtgggaagccagagccgctgcatgccccagccagtcacca
cagcataa