

VTI1B-vesicle transport through interaction with t-SNAREs homolog 1B (yeast) Gene

Reference: PTXBC003142

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VTI1; VTI1-LIKE; VTI1L; VTI2; v-SNARE; vti1-rp1; vesicle transport through interaction with t-SNAREs homolog 1B; vesicle transport v-SNARE protein Vti1-like 1; vesicle-associated soluble NSF attachment protein receptor; vesicle transport through interacti		

For research use only.

GENE DESCRIPTION	vesicle transport through interaction with t-SNAREs homolog 1B (yeast)
GENE ID	10490
ORIGIN SPECIES	Human
ACCESSIONS	BC003142
VECTOR	pDONR223
NCBI SYMBOL	VTI1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctctccgcccctctcggagcatttcgagaagctgcacgagatcttcggcctccatgaagacctacaaggggtgccgagcggtctggggacggcggggacccaagaaaagaagaattgatcagggtatttgataaaagcaacagggaagcaaatgaaacgctg
gcagagatggaggagagctacgttatgcacccctgtcttccgaaaccccatgatgtctaagcttcgaaactaccggaaggaccttgaactccatcgaggaggtgagaagcacacctttgacagccacacctggagccgaggagacatgaatatggcatatctgttagagaatg
agcatatgaatcgctacagtcctcaaaagggaatgctctcaggggcaactgaaagcctgaaccgggcccacccaagattgaacgttctcatcgattgccacagagactgaccagattggctcagaaatcatagaagctgggggaacaacgagaccagttagaacgtaccaagagt
agactggtaaacacaagtgaacttgagcaaaagtcggaagattctcgttcaatgtccagaaaagtgacaaccaacaagctgctgtttccattatcatcttactggagctcgccatcctgggagggcctggttactacaattcttcgcagccattga