

POLR2C-polymerase (RNA) II (DNA directed) polypeptide C, 33kDa Gene

Reference: PTXBC003159

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPB3; RPB31; hRPB33; hsRPB3; DNA-directed RNA polymerase II subunit RPB3; DNA-directed RNA polymerase II 33 kDa polypeptide; DNA-directed RNA polymerase II subunit C; RNA polymerase II subunit 3; RNA polymerase II subunit B3; RPB33; polymerase (RNA) II (D		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide C, 33kDa
GENE ID	5432
ORIGIN SPECIES	Human
ACCESSIONS	BC003159
VECTOR	pDONR223
NCBI SYMBOL	POLR2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtacgccaaccagcctaccgtgcggatcacggagctcactgacgagaatgtcaagttcatcgcgagaacaccgacctggcggtggccaattcgattcggagggtcttcctgctgaggttccataatagccattgactgggttcagattgatccaattcctcagttctcatgatg
aatcattgtctcagagcttggaattatccctcattagtgatgacattgtggacaagctcagttactctcggaactgcacatgtgaggagttctgccccgagtgctcggtggagttcacctcgatgtcggtgcaatgaagaccagacgacatgtcgtcgtgagacctatctccaac
agccccgggtcattccggtgacatccgggaaccgagataatgaccccaatgactacgtggagcaggatgacatctctatcgtcaagttgagaaggccaggagctgagacttcgagcctatgccaaaagggtttggcaaggagcatccaagtgaacctactgcaggggtgg
cttttgaatagatccagacaatgccctgaggcacacagtgtagcccaagcccgaggaatggccaagagtgagtactcggagctggatgaggatgagtcgcaggtccctatgaccccaacggcaagccagaaagtttactacaatgtggagctctgtgctctctgcctgaaac
cattgtcctgtcagccctcaggatgaagaagaactgagtgattacaaactcaattaagccacgagatccagagtgatgtgctaaccataaataa