

GINS2-GINS complex subunit 2 (Psf2 homolog) Gene

Reference: PTXBC003186

| | | | |
|----------------|---|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | HSPC037; PSF2; Pfs2; DNA replication complex GINS protein PSF2; GINS complex subunit 2 (Psf2 homolog); GINS complex subunit 2 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | GINS complex subunit 2 (Psf2 homolog) |
| GENE ID | 51659 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC003186 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | GINS2 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atggacgtccgaggtcgaattcctgccgagaaggagctggtaccattatcccaactcagctctggacaagatctacctatcgggggggactctgggcttttaaccctggttaccctggaagtccccctgtgctggcgattaacctgaacaaagacagaaatgtcgcctgctccctccagatggtgafgafgagaaaagtggagaagatgaggatcatgaacgaaaggagaacttttaccctaatgccagcccttactacatggaacttacgaagctctgttaaatcatgcttcagacaacatcccgaagcagacgaaatccggacctgtcaaggatattgggacactgtatagccaactccgagtgctgctgacagctttgtgagacagcaggaggacatgccaaagctggataactgacctgatggagatcaacaccagcggactttctcacacaagcgtcaaccatgtacaactccgacaaacctcagcctctggagagtactcagctcaggactctag