



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C20orf43-chromosome 20 open reading frame 43 Gene

Reference: PTXBC003359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0549 protein C20orf43; C20orf43; CDAO5; HSPC164; SHUJUN-3; protein RTF2 homolog; replication termination factor 2 domain-containing protein 1; replication termination factor 2 domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 43
GENE ID	51507
ORIGIN SPECIES	Human
ACCESSIONS	BC003359
VECTOR	pDONR223
NCBI SYMBOL	C20orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggttgcagcgggggacaactcccgaaggcgcataactctgtgaagggggccgaagaaggttgagaaggcgcacaagaatgctgaattatggcccaatggaaactattgactctaaagtcaggaaatattaagacaccaaatagttgctgtgaacttgcagactttataacaagaatgcgctcattgaattctctgcagacaactcgcagaagaagctctctgggaaggcgcacatccacatataaagcattaaagaatgtgacagagctgaagctttctgataatctcctctgggaaggggataaaggacaacataaagctgcagacagcagatgacctcagcggggcgctttcattgcccggttggcgctctggagatgacacggcgcagacggctgtctctctctggctcggctgtctgtctctgcagcagcagctgacagaagataaagccgaagttgtccacacgttggcgctctctccagggagatgatgtcatctgtcctaatggcacaaggaggaagtgagcgtctgcacacaagatggagggaggaagagcgtgcagagcgtgcgaagcgtgcgaagcgtgcgaagagcgcgaagacccaagccagagcagctgtgttcaaaacagatgctgtgaagaagcccgagcccatcaaaagttaagacaggggaagcctgcattgtactagagagaagaaaccacgttgcctcccaaaagcacagcaatgaatgagagctctctgcgaagcgtgggaagcctcgtgtgagccacaagaagagctccatctgctgcagctgaagaatcgcgaggcctacaagctccctctttaccactcacagctctcgccaaagcgtcccaaggaggagctgcccaactgggttcacccacactgctcactctctctca

More info: Please contact us for any question: contact@proteogenix.fr