

POMP-proteasome maturation protein Gene

Reference: PTXBC003390

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf12; HSPC014; PNAS-110; UMP1; proteasome maturation protein; 2510048O06Rik; hUMP1; protein UMP1 homolog; voltage-gated K channel beta subunit 4.1; voltage-gated potassium channel beta subunit 4.1		

For research use only.

GENE DESCRIPTION	proteasome maturation protein
GENE ID	51371
ORIGIN SPECIES	Human
ACCESSIONS	BC003390
VECTOR	pDONR223
NCBI SYMBOL	POMP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgccagaggacttgatctgagctaaaggacagtattccagttaacttgcagcaagtgacatttgaagatcatgctcttcggaagggtttctgtgtgaaaaatgaactttgcctagtcacccctgaattatcagaaaaattccagctcaaccaagataaaatgaatttccacactgagaacattcagggtctatttgcctcgtctaaataacagatggaattcaaggcagtcgacgaggttcagcgtcttcattcttcaagctcaaatcttcactggatgtttgaggggtaafgatgagactattggatttgaggatattcttaatgatccatcacaagcgaagtcagggagagccacacttgatggtggaataataaacttggtttactgtaa