



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MOBKL1B-MOB1, Mps One Binder kinase activator-like 1B (yeast) Gene

Reference: PTXBC003398

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MOBKL1B; C2orf6; MATS1; MOB1; MOBK1B; Mob4B; MOB kinase activator 1A; MOB1 Mps One Binder homolog A; MOB1, Mps One Binder kinase activator-like 1B; mob1 alpha; mob1 homolog 1B; mps one binder kinase activator-like 1B		

For research use only.

GENE DESCRIPTION	MOB1, Mps One Binder kinase activator-like 1B (yeast)
GENE ID	55233
ORIGIN SPECIES	Human
ACCESSIONS	BC003398
VECTOR	pDONR223
NCBI SYMBOL	MOBKL1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctctctctcagcagccgctctcttaaacatctaaacaaagaagaataatccctgaagctctcatcagtatgaactcttaaacatcgagagaacactctaggaaagtggaactctgagacaagctgttatgttgcctgaggggagagatctcaatgaatggattctgtgaacactgtggaattcttaaacagatcaacatgttatatggaactattacagaattctgcactgaagcaagctgtccagctcatgtctgcaggctcgagatatgaatactctggcgagatgtactataatataaaaaaccaactcaaatgttctgcaccaaaaactatgactatttgaatgacttgggttcaagatcagcttgaatgaagaaactcttctctctcaagattgggtgtccattctccaaaaactatgagctgtctgagcaagcaactattctaaagctgtctgtcagggttgcctatcttataaacaggaactcatttgaattctgtgagtcagctcgagaagagggcccactcacaacctctttaaagcaacttatttcttctgttcagaaagttaattactctgtatgttgcctgtgacttctcaagaattataaagaagaactctgatacagaagacagata

More info: Please contact us for any question: contact@proteogenix.fr