

**TAF9-TAF9 RNA polymerase II, TATA box binding protein (TBP)-associated factor,
32kDa Gene**

Reference: PTXBC003400

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TAF9 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 32kDa; MGC:5067; STAF31/32; TAF2G; TAFII-31; TAFII-32; TAFII31; TAFII32; TAFIID32; transcription initiation factor TFIID subunit 9; RNA polymerase II TBP-associated factor subunit G		

For research use only.

GENE DESCRIPTION	TAF9 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 32kDa
GENE ID	6880
ORIGIN SPECIES	Human
ACCESSIONS	BC003400
VECTOR	pDONR223
NCBI SYMBOL	TAF9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctctggcaagacgcttctcccaagagcgcgaagatgcacagatgatggcacaatcctgaaggatattgggattacagaatgatgccaagagtataaatcagatgttgagtttgcctccgatatgtgaccacaattctagatgatcaaaaattattcaagccatgctaa
gaaagctactgttgatgcagatgatgtgcgattggcaatccagtgccgcgctgatcagctctttacctctctcccccaagagatttttattagatattgcaaggcaagaatcaaaccccttgcattgatcaagccatattcaggctcctaggttgcacactgatagatactgcttaacagctcc
aaactataggctgaaatctttacagaaaaaggcatcaacttctgcgggaagaataacagtcgccggttaagtgtgttcagttactagcagaccaagtaactccacactaggcacaccaacccacagaccatgtctgtttcaactaaagtagggactcccatgtccctcacaggtcaaa
gtttacagtcagatgcctactctcagctcctcagctgttaaaagcttcaattcctgcaacctcagcagttcagaatgttctgattaatccatcattatcggttccaaaacatttttaccactaatatgatgtcatcaaaaatactgccaatgaatcatcaaatgcattgaaaagaaacgtgaa
gatgatgatgatgacgatgatgatgatgactatgataatctgtaa