

RAP2C-RAP2C, member of RAS oncogene family Gene

Reference: PTXBC003403

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAP2C, member of RAS oncogene family; small GTPase RAP2C; ras-related protein Rap-2c		

For research use only.

GENE DESCRIPTION	RAP2C, member of RAS oncogene family
GENE ID	57826
ORIGIN SPECIES	Human
ACCESSIONS	BC003403
VECTOR	pDONR223
NCBI SYMBOL	RAP2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggaatacaagtagtggttagggagtgagggtggcgaatctgccttactgtgcagttgtcactgggactttcattgagaaatgatgacccaccattgaagatttaccgcaagagatcgaagtggactttccccctcgtgctggaaattctggacaccgcaggaactg
agcagtttgccctcatgagagatctctacatcaaaaacggccaaggtttcatcctggtttatagcctggttaatacaacagctctttcaggatatcaagccaatgagagatcaaatgtcagagtgaagagatatgaaaagtcccactaatcctagtaggaaataaagtggatctggaaccagaaa
gagaggttatgtcttcagaaggcagagctctggttcaagaatggggctgtcctttcatggagacatcggcaaaaagtaaatcaatggtggaactttttctgagatcgtcaggcaaatgaactattcatccctgccggagaagcaagatcagttgttacaactgtgtcgtccagtaa