

CXCL14-chemokine (C-X-C motif) ligand 14 Gene

Reference: PTXBC003513

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMAC; KEC; KS1; MIP-2g; MIP2G; NJAC; SCYB14; C-X-C motif chemokine 14; CXC chemokine in breast and kidney; MIP-2 gamma; bolekin; breast and kidney; chemokine (C-X-C motif) ligand 14; chemokine BRAK; small inducible cytokine subfamily B (Cys-X-Cys), membre		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 14
GENE ID	9547
ORIGIN SPECIES	Human
ACCESSIONS	BC003513
VECTOR	pDONR223
NCBI SYMBOL	CXCL14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccctgctccacgccgcgccctcgggtcagcatgaggtcctggtggcggcgctgctgctgctggtgtacaccgcggtgtggacgggtccaaatgcaagtgtctccggaaggggaccaccaatccgtacagcgacgtgaagaagctggaatgaagccaaag
taccgcactgcgaggagaagatggttatcatcaccaccaagagcgtgtccaggtaccgaggtcaggagcactgcctgcacccaagctgcagagcaccgaagcgttcatcaagtgtgtacaacgcctggaacgagaagcgcgagggtctacgaagaatag