

### WDR32-WD repeat domain 32 Gene

Reference: PTXBC003520

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WDR32; DDB1- and CUL4-associated factor 10; WD repeat domain 32; WD repeat-containing protein 32; DDB1 and CUL4 associated factor 10		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	WD repeat domain 32
<b>GENE ID</b>	79269
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC003520
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	WDR32
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

#### ORF (Only for reference, product provided without STOP codon)

atgcgaatgaggttaacaccagattgtccaaaatgttgaattcaacgtcctctgatatctttaaatttgcacgtgacttaactaactctttagaagtagcagctatccatttaagagcaagaagaactactcaagtcaggagttcaccacgaaatagcttgaagtcgtaaccctga  
agttcttggfagagtgaccacggaaactgcatcacgtccttacagctgcacccaaaaggctggccactcttctcggfctcaagcaacagtgatgatgagagtgacttggctatgaattccaagaaggagctccagtcgacgtgtctcactctcgtttctctacgactgactcatta  
cattgaagaagccaatgtaggcagggttacatcaagaactttgcttcagccccgatggacgaatgattctccccacatggctatgggattcgttgggattgacaacagtcagtgaaactgttgactccttccccaaagaagccagtcctcctgaggatcctctctgtactc  
tcataatgatgtgtactgacaacaaagtctccaacacatgtcagattgcctcagggtgccttagtgacgggtttcttctatcagccaaagtttag