



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

POLL-polymerase (DNA directed), lambda Gene

Reference: PTXBC003548

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BETAN; POLKAPPA; DNA polymerase beta-2; DNA polymerase beta-N; DNA polymerase kappa; polymerase (DNA directed), lambda; polymerase (DNA) lambda		

For research use only.

GENE DESCRIPTION	polymerase (DNA directed), lambda
GENE ID	27343
ORIGIN SPECIES	Human
ACCESSIONS	BC003548
VECTOR	pDONR223
NCBI SYMBOL	POLL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatcatcaaaaagtactctgcaaaagttctcaggaagggagagggagaaagcagaaggaggcctgcgatctccctgggattgggaagcggatgctcgagaaaatcatagagatcctggagagcgggcatttgcggaagctggaccatcatcagtgagagcgtgctgtctggagc
tcttctccaaactctggggagctggcgacgaagatgccagatgtgtgtccaaacagcggcttgcgatctcggaagacatcgcgacgagcctccttcgacaacccacagcggcattccgaagcattacagtgacttcttggacgtgtccggagggagggatcgacacagattga
cgacagatgcacgaagaacgcggcagcggcttctaactctgggctctgtgtgtggcatctgtctacacgcgggaagagcgactctgtgtgatgtgcagctgtctcatcaccacacagatggcgcgttcaccgggttcgcgcctctctgacagtctgcgcgaaggggttc
ctcacagatgactgtgtgacccaagaggaatggtccacacagaagtactctgggggtgtgccggctccacggccagggcggcgagcgcgcctggacatcatcgtgtgtgccatcagcagatttgggtgtgcctctacttaccggctctgcacacttcaaccgcctcatg
cagccctggcccaaaccaaggcgcgatgctgtcgaacatgcctctcagcactgtgtgtgtccggaacaccatcggctgcgaagctggggcctgcccgaagtgtctgccctccactgagaaggatgtcttcaggctcttagcgctccctaccgcgaagacgtctgacgcgggactgtg
ga

More info: Please contact us for any question: contact@proteogenix.fr