



Reference: PTXBC003551

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TG(C); TGC; protein-glutamine gamma-glutamyltransferase 2; C polypeptide, protein-glutamine-gamma-glutamyltransferase; TGase C; TGase H; TGase-2; tissue transglutaminase; transglutaminase C; transglutaminase H; transglutaminase 2		

For research use only.

GENE DESCRIPTION	transglutaminase 2 (C polypeptide, protein-glutamine-gamma-glutamyltransferase)
GENE ID	7052
ORIGIN SPECIES	Human
ACCESSIONS	BC003551
VECTOR	pDONR223
NCBI SYMBOL	TGM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr