



Reference: PTXBC003560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RIBIIR; RPN-II; RPNII; SWP1; dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit 2; dolichyl-diphosphooligosaccharide--protein glycosyltransferase 63 kDa subunit; oligosaccharyltransferase complex subunit (non-catalytic); ribophorin-2;		

For research use only.

GENE DESCRIPTION	ribophorin II
GENE ID	6185
ORIGIN SPECIES	Human
ACCESSIONS	BC003560
VECTOR	pDONR223
NCBI SYMBOL	RPN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr