

TCL1A-T-cell leukemia/lymphoma 1A Gene

Reference: PTXBC003574

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T-cell leukemia/lymphoma protein 1A; T-cell lymphoma-1; oncogene TCL-1; oncogene TCL1; protein p14 TCL1; T-cell leukemia/lymphoma 1A		

For research use only.

GENE DESCRIPTION	T-cell leukemia/lymphoma 1A
GENE ID	8115
ORIGIN SPECIES	Human
ACCESSIONS	BC003574
VECTOR	pDONR223
NCBI SYMBOL	TCL1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgagtgcccgacactcggggaggcagtcaccgaccacccggaccgctgtggccctgggagaagtcgtgtatttgacgagaagcagcagcctggctgcccttaaccatcgagataaaggataggttacagttacgggtgctcttgcgtcgggaagacgtcgtcctgggga
ggcctatgacccccaccagataggcccaagcctgtgcctatcatgtggcagctctaccctgatggacgataccgatcctcagactccagtttctggcgttagtgtaccacatcaagattgacggcgtggaggacatgcttctcgagctgctccagatgactga