

URM1-ubiquitin related modifier 1 homolog (S. cerevisiae) Gene

Reference: PTXBC003581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf74; ubiquitin-related modifier 1; ubiquitin-related modifier 1 homolog; ubiquitin related modifier 1		

For research use only.

GENE DESCRIPTION	ubiquitin related modifier 1 homolog (S. cerevisiae)
GENE ID	81605
ORIGIN SPECIES	Human
ACCESSIONS	BC003581
VECTOR	pDONR223
NCBI SYMBOL	URM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcgccctgtcagtgaggtgagttcggaggtggtgcggagctcctgtttgacggtattagaacacatcagtcactttgcctggacaggaggaaacctgggacatccggaacctgctcatctggatcaagaagaatttgctaaaagagcggccagagttgtcatccaggag
acagcgtgcggccagggaattctgtgctgattaacgatgccgactgggagctactgggtgagctggactaccagcttcaggaccagacagcgtcctcttcctccactctgcacggcggtga