

C21orf33-chromosome 21 open reading frame 33 Gene

Reference: PTXBC003587

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ES1; GT335; HES1; KNPH; KNPI; ES1 protein homolog, mitochondrial; Keio novel protein I; human HES1 protein, homolog to E.coli and zebrafish ES1 protein; testis secretory sperm-binding protein Li 237E; chromosome 21 open reading frame 33		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 33
GENE ID	8209
ORIGIN SPECIES	Human
ACCESSIONS	BC003587
VECTOR	pDONR223
NCBI SYMBOL	C21orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctgtgagggccctggtgacctgaggtcgtgcggcatctgcattcacgtccctgtccccggcggtcggagcctccagcgcgcagccctcacctctcgtgccgcgcccgccgaggggtcgcgtggtgtgtctggtgatgcggagtcacgatgggaccgaga
tcacagaggcctcggcgatcctggtgcacctgagccgtggagggtgtaagtccagatcttgcctctgacgtccctcagatgcacgtgattgaccacacaaaggcgagccgtccgaaggcgagagcaggaatgtttgaccgagtcgcgaggtatgccctgtggcaaatcacaga
cctggccaacctcagtcagccaacctgatgctgccatcttccaggagccttggagcggcctaaaaacctgagcacgtttgccgtggacgggaaagattgcaaggtgaataaagaagtgagcgtgtcctgaaggagttccaccaggccgggaagcccatcggtgtgtcattgc
acctgtctcgcggccaaggtgctcagaggcgtcaggtgactgtggccacgagcaggaaggtggcaagtgccctatgccgggaccgcagagccatcaaggccctgggtgccagcactgcgtgaaggaggtgtcgaagctcacgtggaccagaaaaaacagggtgt
cacgacccagccttcattgtgcgagacggcactccactacatccatgatggatcggagccatggtgaggaaggtgctggaactacgtgaaagtga