

MRPS25-mitochondrial ribosomal protein S25 Gene

Reference: PTXBC003590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP-S25; RPMS25; 28S ribosomal protein S25, mitochondrial; S25mt; mitochondrial 28S ribosomal protein S25; mitochondrial ribosomal protein S25		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S25
GENE ID	64432
ORIGIN SPECIES	Human
ACCESSIONS	BC003590
VECTOR	pDONR223
NCBI SYMBOL	MRPS25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatgaaggcgcttccccatccgcaccctgcaatatctgagccagggaacgtggtgttcaaggactccgtgaaggtcatgacagtgaattacaacacgcatggggagctggcgaggcgccaggaagttgtgttttcaacatacctcagattcaataaaaaacctt
gggtgcagatcatgatgtttaagaacatgacgcccgtcaccctctcgtgattctacttagattctggggagcaggtctggtggtgagacgaagcaataaggagatcatggagcacatcagaaaaatcttggggaagaatgaggaacccctcagggaaggaggagagaaa
aagcagctttctaccagccaacttcggccctgaaagtactgcctgcgggagtgcatctgtgaagtgaaggcaggtgcccctgccccagcctgggtgccattaccaaggagatgaggggaagtacaaagccgtctgaaagccgatgccagactaa